

System Requirements of Scilligence Products

Most our products can be deployed on the cloud or/and in-premise environment.

The system parameters specified below are recommended for successful product deployment. Our products can be deployed on SQL Server Cluster and IIS Web for medium to large size organizations.

Product	Windows Server	CPU	Memory	IIS/ASP.NET	SQL Server	or Oracle	SharePoint
Scilligence ELN	-2008 R2 x64 -2012 x64	1+ Quad Core	6+GB	.NET 2.0	-2008 R2 x64 -2012 x64	9, 10 or 11	
Chem4SharePoint	-2008 R2 x64	1+ Quad Core	8+GB	.NET 3.5	-2008 R2 x64	9, 10 or 11	2010
Chrawler	-2008 R2 x64 -2012 x64	1+ Quad Core	6+GB	.NET 3.5	-2008 R2 x64 -2012 x64	9, 10 or 11	
RegMol	-2008 R2 x64 -2012 x64	1+ Quad Core	4+GB	.NET 2.0	-2008 R2 x64 -2012 x64	-	
MolSQL	-2008 R2 x64 -2012 x64	1+ Quad Core	4+GB	.NET 2.0	-2008 R2 x64 -2012 x64	-	
MolEngine	-	-	-	.NET 2.0	-	-	
JSDraw.WebServices MolEngine.WebServices	-2008 R2 x64 -2012 x64	1+ Quad Core	4+GB	.NET 2.0	-2008 R2 x64 -2012 x64	-	

Last Modified on January 27, 2013